

VERSAWAVE™ ACCESSORIES

Choose the one you need for your practice

Unlike the drill, there is little tactile feedback with laser dentistry since a beam of light is doing the work. The tips offered by HOYA ConBio maximize your ability to work quickly and precisely. Maximum visibility is provided by a patented, thin curved design which emits the water spray right at the target, eliminating overspray. This is especially helpful when working with a mirror on the upper arch.

A right angle handpiece sleeve is also offered when access in the posterior is more important.

These accessories are combined with the Perfect Balance™ Fiber Support System to provide you with a zero-resistance handpiece that you can use without effort all day long.

Features:

- Unique patented curved tips
- Quartz and sapphire tips
- Unique water spray at tip end
- Flexible endodontic tips
- 360° swivel, lightweight and compact handpiece
- All accessories autoclavable and reusable

To Order Call

Lions Dental Supply
951-276-3225

HOYA ConBio
HOYA & CONBIO

©2005 HOYA ConBio

► The Master's Series of Hard and Soft Tissue Tips



STRAIGHT TIP
600 MICRON

30°
600 MICRON

80°
600 MICRON



SOFT TISSUE
CONICAL

80°
400 MICRON

PERIO 90°
400 MICRON

► Prep-Pro™ Sapphire Tips



600 MICRON

700 MICRON



TAPER
450 MICRON

CHISEL
1200 x 400 MICRON

► The Endodontic Series



200 MICRON
31 MM

300 MICRON
31 MM

► Exclusive Handpiece Sleeves



STRAIGHT SLEEVE

RIGHT ANGLE

SUPERB ENGINEERING LEGENDARY RELIABILITY