

## Cranberry SIGMA 5810 Series

Cranberry SIGMA Powder-Free Latex Examination Gloves



- New Improved Formula
- Powder-Free
- Available Sizes: XS,S,M,L, & XL
- \*Pro2S Matrix: **Powder Free, Barrier Protection**

### Solution Profile

#### *Solution Type*

- Material: **Quality Natural Rubber Latex**
- Fitting: **Ambidextrous**
- Cuff: **Beaded**
- Surface: **Textured**

#### *Ideal Solution for*

- Odorless Preference
- Powder Irritation
- Enhanced Grip
- Good Tear Resistance

**Solution Highlights**

- Excellent grip for usage in wet and dry situations.
- Beaded cuffs for better tear resistance.

**Quality Standards**

- **Surpasses United States ASTM D3578-91 standard for medical grade examination gloves**
- **Complies with American National Standard/ American Dental Association (ADA) Spec. No. 76-1991 ISO 9002**
- **Certified by Standards and Industrial Research Institute of Malaysia (SIRIM).**
- **Endorsed by TUV, Germany**

**Ordering Numbers**

- 5815 (Extra Small)
- 5816 (Small)
- 5817 (Medium)
- 5818 (Large)
- 5819 (Extra Large)

**Technical Specifications**

|  | ASTM***  | SIGMA     |
|--|----------|-----------|
| Length   | min. 230 | min. 240  |
| Finger Thickness (mils)                            | 3.15     | min. 0.15 |
| Palm Thickness (mils)                              | 3.15     | min. 0.12 |
| Tensile Strength (MPa)                             |          |           |
| Before Aging                                       | min 21   | 14        |
| After Aging  | min 16   | 14        |
| Elongation (%)                                     |          |           |
| Before Aging                                       | min. 700 | min. 500  |
| After Aging  | min. 500 | min. 700  |
| ***ASTM D 5250. Test Result From Production Sample |          |           |

**Ecological Information**

- Bio-Compatibility: The chemical formulation of the gloves and surface lubricating materials do not contain any substances normally known to be harmful to the user or to any person with whom the gloves contact.

**Other Information**

- To prolong product life, gloves should be stored in a cool, dry location, located away from sources of heat and light.

**Packaging**

- 100 gloves/box by weight(90 for XL)
- 10 boxes/case