

Cranberry SILKCARE 7810 Series Cranberry SILKCARE Powder Free Latex Examination Gloves



- Microweb Textured
- 50ug/g water extractable protein
- Protein-Rated to Minimize Allergen Content , FIRST TOUCH™ Enhanced
- Available Sizes: XS,S,M,L, & XL
- *Pro2S Matrix: **Powder Free, Barrier Protection**

Solution Profile

The Data in this Material Safety Data Sheet relates to the specific material designated herein and does not relate to the use in combination with other material or in any process.

Solution Type

- Material: **Quality Natural Rubber Latex**
- Fitting: **Ambidextrous**
- Cuff: **Beaded**
- Surface: **Smooth**

Ideal Solution for

- Effective Barrier Protection
- Improved Donning Efficiency
- Good Puncture Resistance

Solution Highlights

- Moisturized inner coating offers superior donning and soothing comfort.
- Exceed worldwide standards for barrier protection, strength, comfort and low protein to minimize the risk of latex protein allergy.
- Specially coated with USP grade Lanolin and Vitamin E, designed to minimize dry skin irritation associated with frequent hand washing.
- The moisturizing ingredients are formulated to soften without leaving an oily residue.

Quality Standards

- **United States ASTM D3578 standard for medical examination gloves**
- **American National Standard/American Dental Association (ADA) Spec. No. 76**
- **TUV Certification, Germany**
- **ISO 9002 Certified**

Ordering Numbers

- 7815 (Extra Small)
- 7816 (Small)
- 7817 (Medium)
- 7818 (Large)
- 7819 (Extra Large)

Technical Specifications

	ASTM***	Silkcare
Length	min. 230	min. 240
Finger Thickness (mils)	3.15	min 6.5
Palm Thickness (mils)	3.15	min 6
Tensile Strength (MPa)		
Before Aging	min 21	24
After Aging	min 16	20
Elongation (%)		
Before Aging	min. 700	min. 800
After Aging	min. 500	min. 800
***ASTM D 3578.SMG 1007.ISO9002.European CE.Test Result From Production Sample		

Ecological Information

- Bio-Compatibility: The chemical formulation of the gloves and surface lubricating materials do not contain any substances normally known to be harmful to the user or to any person with whom the gloves contact.

Other Information

- To prolong product life, gloves should be stored in a cool, dry location, located away from sources of heat and light.

Packaging

- 100 gloves/box by weight(90 for XL)
- 10 boxes/case
- Protein Levels
This glove contains 50 micrograms or less of total water extractable protein per square decimeter of glove surface area. Total protein content is determined by using Modified Lowry Test Method, ASTM D5712-99.
- Latex Caution Statement
Caution: This product contains natural rubber latex which may caused allergic reactions.