


Cranberry **Xlim™** Series  
Cranberry Xlim Powder Free Nitrile Examination Gloves



- Ultra Soft Nitrile
- 100% Latex-Free
- Odorless
- "First Touch" Enhanced
- Available Sizes: XS,S,M,L, & XL
- \*Pro2S Matrix: **Latex Free, Powder Free, Barrier Protection**

### Solution Profile

#### *Solution Type*

- Material: **Nitrile**
- Fitting: **Ambidextrous**
- Cuff: **Beaded**
- Surface: **Textured**

#### *Ideal Solution for*

- Latex Allergy
- Powder Irritation
- Good Puncture Resistance
- Enhanced Grip and Lotion Compatibility

#### *Solution Highlights*

- **New proprietary slim formulation provides "2nd skin" comfort and delivers unbelievable tactile sensitivity.**

- Compliment our existing range of Contour and Contour Plus nitrile gloves.
- Available in exclusive cool blue color and compact packaging for easy storage.
- Competitively priced for customers looking for a comfortable synthetic solution.

**Quality Standards**

- **United States ASTM D3578-00 standard for medical examination gloves**

**Ordering Numbers**

- 3125 (Extra Small)
- 3126 (Small)
- 3127 (Medium)
- 3128 (Large)
- 3129 (Extra Large)

**Technical Specifications**

	ASTM***	Contour
Length	min. 230	min. 235
Finger Thickness (mils)	3.2	min 5.5
Palm Thickness (mils)	2.0	min 4.5
Tensile Strength (MPa)		
Before Aging	min 9	17
After Aging	min 9	15
Elongation (%)		
Before Aging	min. 300	min. 700
After Aging	min. 300	min. 500
***ASTM D 3578-00, 5712-99, 6124-00. EN 455 part 1,2. FDA 1000ml Watertight Test. Test Result From Production Sample		

**Packaging**

- 100 gloves/box by weight
- 10 boxes/case