

clearly nitrile™

Cranberry Series
Cranberry Clearly Nitrile Powder Free Nitrile Examination Gloves



- Ultra-Soft Comfort
- 100% Latex-Free
- Odorless & Colorless
- SmartGrip™ Enhanced
- Available Sizes: XS,S,M,L, & XL
- *Pro2S Matrix: **Latex Free, Powder Free, Barrier Protection**

Solution Profile

Solution Type

- Material: **Nitrile**
- Fitting: **Ambidextrous**
- Cuff: **Beaded**
- Surface: **Smooth**

Ideal Solution for

- Latex Allergy
- Powder Irritation
- Good Puncture Resistance
- SmartGrip™ for dry or wet donning

Solution Highlights

- The only even surface nitrile glove with SmartGrip™ Technology designed to deliver unparalleled tactile sensitivity and optimal grip (dry or wet)
- Compliments our existing range of Contour Plus, Contour, and Xlim nitrile gloves.
- Odor Free, Color Free and Latex Free.
- Ultra soft comfort with superior protection.

Quality Standards

- **United States ASTM D3578-00 standard for medical examination gloves**

Ordering Numbers

- 3325 (Extra Small)
- 3326 (Small)
- 3327 (Medium)
- 3328 (Large)
- 3329 (Extra Large)

Technical Specifications

	ASTM***	ClearlyNitrile™
Length	min. 230	min. 240
Finger Thickness (mils)	3.2	min 6.0
Palm Thickness (mils)	2.0	min 5.0
Tensile Strength (MPa)		
Before Aging	min 9	14
After Aging	min 9	14
Elongation (%)		
Before Aging	min. 300	min. 500
After Aging	min. 300	min. 400
***ASTM D 3578-00, 5712-99, 6124-00. EN 455 part 1,2. FDA 1000ml Watertight Test. Test Result From Production Sample		

Packaging

- 100 gloves/box by weight
- 10 boxes/case