



## POWDER-FREE GREEN LATEX EXAMINATION GLOVES



### GLOVE FEATURES

Length (mm)	240 min
Cuff Thickness (mm)	0.11 ± 0.03
Palm Thickness (mm)	0.14 ± 0.03
Finger Thickness (mm)	0.15 ± 0.03
Tensile Strength	
Before Aging (MPa)	21 min
After Aging (MPa)	16 min
Elongation	
Before Aging (%)	700 min
After Aging (%)	500 min
Protein Rating	max 50µg/g
Powder Residue	max 2.0 mg/glove
Quality Assurance Conformance	ASTM 357800 SMG 1007 ISO 9001 European CE

### GLOVE DESCRIPTION

Material	Natural Rubber Latex
Color	Green
Scent	Citrus Mint
Cuff	Beaded
Interior Surface	Smooth polymer interior coating Powder-free
Exterior Surface	Micro-textured surface
Usage	Single use only Fits either hand
Size	X-Small, Small, Medium, Large
Packaging Configuration	100 gloves by weight/box 10 boxes/carton

### ORDERING INFORMATION

Item Number	Size
2835	X-Small
2836	Small
2837	Medium
2838	Large

- Refreshing citrus mint scent & winter green color provides pleasant wearing experience.
- Smooth polymer interior coating improves moist or dry donning.
- Micro-textured surface deliver consistent grip under all operating conditions.

#### **Caution Statement:**

Safe use of this gloves by or on latex sensitized individuals has not been established.



**Lions Dental Supply**

**951-276-3225**

**[www.lionsDentalSupply.com](http://www.lionsDentalSupply.com)**

