



POWDER-FREE LATEX EXAMINATION GLOVES



GLOVE FEATURES

Length (mm)	240 min
Cuff Thickness (mm)	0.08
Palm Thickness (mm)	0.10
Finger Thickness (mm)	0.12
Tensile Strength	
Before Aging (MPa)	18 min
After Aging (MPa)	14 min
Elongation	
Before Aging (%)	650 min
After Aging (%)	500 min
Protein Rating	max 50µg/dm ²
Powder Residue	max 2.0 mg/glove
Quality Assurance Conformance	ASTM D 3578 Latest revision ASTM D 5712 Latest revision ASTM D 6124 Latest revision ISO 11193 EN 455 Parts 1 & 2 FDA 1000ml Watertight Test

- Smooth polymer interior coating improves moist or dry donning.
- Micro-textured surface deliver consistent grip under all operating conditions.

GLOVE DESCRIPTION

Material	Natural Rubber Latex
Color	Natural
Cuff	Beaded
Interior Surface	Smooth polymer interior coating Powder-free
Exterior Surface	Micro-textured surface
Usage	Single use only Fits either hand
Size	X-Small, Small, Medium, Large
Packaging Configuration	100 gloves by weight/box 10 boxes/carton

Caution Statement:

Safe use of this gloves by or on latex sensitized individuals has not been established.



ORDERING INFORMATION

Item Number	Size
2815	X-Small
2816	Small
2817	Medium
2818	Large

Lions Dental Supply

951-276-3225

www.LionsDentalSupply.com

