

# ART.

## ART-E1

# ELECTROSURGERY UNIT

Effective & Reliable



BONART MEDICAL TECH. INC.

# ART-E1 ELECTROSURGERY UNIT

User friendly control panel  
Safely designed electronics  
Minimal scar tissue formation  
Minor to high demand dental surgery

CUT, COG 1 & COG 2 mode  
10 level power setting  
7 pc. Various standard electrodes

## Specifications

Power supply	115 V $\pm$ 5% ~ 50/60 Hz 92 VA 230 V $\pm$ 5% ~ 50/60 Hz 92 VA
Max. output power	50 W ( @ 400 ohm load)
Work frequency	1.5 ~ 1.7 MHz $\pm$ 5%
Dimension	23 (L) x 22.5 (W) x 8.5 (H) cm
Weight	4 Kg
Handpiece cable	190 cm
Ind. plate cable	190 cm

## Included in delivery

- Handpiece x 1
- Ind. Plate x 1
- Electrodes set x 1
- Foot switch x 1
- AC cord x 1
- User' s manual x 1



• CUT, COG 1  
& COG 2 MODE



• ELECTRODES SET



CE ISO-13485  
0434 US 510(k)



**Lions Dental Supply**  
**951-276-3225**  
[www.LionsDentalSupply.com](http://www.LionsDentalSupply.com)

