

ENDOMAX™

CORDLESS ENDODONTIC HANDPIECE



ISO 9001
ISO 13485

GMP

ENDOMAX™


- Cordless convenience offers complete portability.
- Auto-stop and auto-reverse mechanisms for patient safety.
- Heavy-duty high torque motor.
- Reciprocating setting available without use of additional attachment.
- Ergonomic design provides comfort and easy operation.
- Features five torque settings.
- Automatic forward/reverse functions.
- Nine rotation speeds available for different applications.
- High capacity lithium-ion battery provides stable and prolonged usage.
- Offers quick and easy coupling.
- Low battery warning light.

DENTAMERICA®

ENDOMAX™

CORDLESS ENDODONTIC HANDPIECE



Auto-Reverse Forward:



Reciprocating Mode:



Auto-Reverse Stop:



SPECIFICATIONS:

INPUT	AC 100V~240V 50/60Hz
FUSE	500mA 250V
BATTERY TYPE	3.7V 1620mAh
DIMENSIONS (HANDPIECE)	L29xW25xH198mm
DIMENSIONS (CHARGING STATION)	L105xW85xH59mm
WEIGHT (HANDPIECE)	155g
WEIGHT (CHARGING STATION)	125g

Lions Dental Supply

951-276-3225

www.LionsDentalSupply.com

DENTAMERICA®