


The use of CMOS technology in the QuickRay USB sensor allows to make the best use of the geometry of the sensor to maximize the Image Area.

Comfort of use is therefore greater and optimal Image Quality is achieved every time an X-ray is shot

The composition of a QuickRay Direct USB Sensor



But the real difference is the amazingly low cost of ownership for the QuickRay sensor. Not only is QuickRay the most aggressively priced sensor on the market, but it also comes with the most comprehensive and most affordable Warranty plan: The first 2 years of warranty are provided for free and VDC proposes an extended warranty option for 3 more years.

What's more, if you choose not to select VDC's extended warranty, QuickRay Direct USB features a cable that can be replaced for a low fee if broken, pinched or otherwise

Preferred Imaging Software: CHOICE IMAGING

With over 11,000 licenses sold worldwide CHOICE Imaging is a leading Dental Imaging Software.

Video dental Concepts introduces a real breakthrough in Digital X-Ray Sensors with its new line of QuickRay Direct USB sensors. Available in Size #1 and Size #2, QuickRay sensors deliver the Highest Image Quality with its CMOS/ CsI technology providing more than 26 lp/ mm in 4096 grey scale.

The QuickRay rounded corners make positioning easier and more comfortable for the patient. And it's a direct USB sensor: no need for an interface box... directly plug it in the computer and you are ready to take your X-Ray !



www.videodental.com